



INTRODUCTION

- To obtain virus isolates for phenotypic analysis, clinical specimens are generally cultured and titrated in absence of inhibitors prior to further characterization.
- Proportions of virus variants may change during culture, and minority species with reduced sensitivity may remain undetected.
- Here we developed sensitive detection technology for direct phenotypic characterization of virus in samples containing between 1E2 and 1E7 plaque forming unit (PFU)/mL:
 - Within 24h after inoculation, without prior virus culture or titration
 - Facilitating isolation of clonal progeny virus populations for further phenotypic and genotypic analyses.
 - Automated imaging ensures objective data and traceability of infected cell patterns which formed in presence and absence of antiviral substances.
 - Format compatible with large scale sample testing


METHODS



RESULTS




- ONE spot showing ONE infected cell
- TOO EARLY to detect secondary infections




- Monolayer with ONE spot showing that progeny virus of ONE infectious unit infected >>1000 cells within 24 h

Sensitive A(H3N2) virus (21h)




Resistant A(H3N2) virus (21h)





% INHIBITION of maximal virus signal (row H)

| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|------|------|-----|-----|-----|-----|-----|------|------|------|------|------|
| A | 93 | 98 | 99 | 99 | 100 | 100 | - | - | - | - | - | - |
| B | 81 | 90 | 98 | 100 | 100 | 100 | - | - | - | - | - | - |
| C | 51 | 76 | 89 | 100 | 100 | 100 | - | - | - | - | - | - |
| D | 46 | 48 | 63 | 90 | 97 | 100 | - | - | - | - | - | - |
| E | 26 | 54 | 43 | 8 | 30 | 65 | - | - | - | - | - | - |
| F | 11 | 42 | 30 | 29 | 5 | 9 | - | - | - | - | - | - |
| G | 5 | 18 | 25 | 41 | -9 | -16 | - | - | - | - | - | - |
| H | 11.5 | 10.4 | 7.7 | 6.7 | 5.4 | 4.7 | 1.8 | 0.00 | 0.10 | 0.00 | 0.00 | 0.00 |



| | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | 11 | 12 |
|---|-----|-----|-----|-----|-----|-----|-----|------|------|------|------|------|
| A | 7 | -1 | -56 | 18 | 42 | - | - | - | - | - | - | - |
| B | 7 | 27 | -62 | -11 | 71 | - | - | - | - | - | - | - |
| C | 31 | 34 | -45 | 14 | 78 | - | - | - | - | - | - | - |
| D | 14 | 35 | -20 | 9 | 33 | - | - | - | - | - | - | - |
| E | -1 | 34 | -25 | -6 | 19 | - | - | - | - | - | - | - |
| F | 8 | 27 | -12 | 0 | 45 | - | - | - | - | - | - | - |
| G | -11 | 30 | -48 | -11 | -19 | - | - | - | - | - | - | - |
| H | 5.0 | 5.0 | 3.2 | 4.5 | 3.9 | 2.2 | 0.1 | 0.01 | 0.00 | 0.01 | 0.00 | 0.00 |



CONCLUSIONS

- This phenotypic ViroSpot™ Assay allows RAPID DETECTION and characterization of virus in samples containing 1E2 – 1E7 PFU/mL.
- Progeny virus of RESISTANT virus clones, which were cultured in the PRESENCE of inhibitor, can be ISOLATED and stored for further analyses.
- Since prior culture in absence of inhibitors is not required, this assay may provide valuable PHENOTYPIC information on HETEROGENEOUS VIRUS POPULATIONS in clinical specimens.